# **Buccal Mucosal Bleeding Time Procedure**

### **Summary**

The Buccal Mucosal Bleeding Time (BMBT) can provide the clinician with important information regarding the interaction between platelets and the vascular endothelium. Although the BMBT can be prolonged in conjunction with a variety of hemostatic disorders (Thrombocytopenia, vascular injury, or hemophilia), a prolonged bleeding time in the presence of a normal platelet count suggests defective platelet function. Platelet dysfunction is commonly seen with von Willebrand's Disease (VWD); however, aberrant platelet function may also accompany other diseases or the administration of certain drugs. A BMBT test should be considered for all patients undergoing surgery, but especially for those breeds known to be prone to VWD.

#### **Materials Provided**

Tenderfoot® Lancet Device\*
Blotting Paper (#1 Whatman is recommended.)
Stopwatch

\*Important: Due to discontinuation of the Surgicutt® line of products, the Tenderfoot® Lancet, a device used for collecting blood samples from young children, has been selected as a suitable alternative. Both Surgicutt® and Tenderfoot® lancets produce a 1 mm depth incision however the length of the Tenderfoot® incision is 2.5 mm vs. 5 mm for Surgicutt®. This should not affect the expected results especially in the case of coagulation-impaired animals. However, when transitioning to the new device, it is strongly recommended that users test 5-10 normal patients to assess any changes to the baselines already established for your practice.

## **Procedure**

- 1. Position the patient in lateral recumbency.
- 2. Firmly tie a gauze muzzle around the proximal maxilla while everting the upper lip to expose the buccal surface.
- 3. Remove the safety guard and place the Tenderfoot® device lightly against the exposed buccal surface and release the spring-loaded lancet. The blade will produce a 1 mm deep linear incision.
- 4. Start the stopwatch and begin absorbing blood away from the incision at approximately 5 second intervals. Blood should be "wicked away" from the edge of the incision. Do not apply direct pressure.
- 5. Continue blotting until bleeding stops. Record the elapsed time.







Remove safety guard



Perform Incision



Blot every 5 seconds



Record elapsed time

### Normal Values\*

Dogs < 4 minutes Cats < 3 minutes

\*Values extracted from the literature. Each clinic should establish its own reference value by testing 5-10 normal samples.

#### **Additional Notes**

Patients with prolonged bleeding times should receive further investigation. Prolonged bleeding times can occur secondary to Thrombocytopenia. A platelet count should accompany the BMBT.

Reorders:

HEM-BMBTKIT Complete Buccal Mucosal Bleeding Time Kit (Surgicutt, Blotting Paper, Stopwatch)

HEM-ACCTF501 Tenderfoot® Lancet; Pkg. 10
GEN-WH1002-110 Whatman Filter Paper (for blotting)

Distributed by: Vetlab Supply · Palmetto Bay, FL 33157 · 800.330.1522 ·Fax 305.232.8421 ·orders@vetlab.com