
## Preparation of Standardized Suspensions for Antibiotic Susceptibility Testing

(McFarland Turbidity Method)


- 1. Prepare test suspension by inoculating a broth tube with 1-3 like colonies harvested from culture plate.
- 2. Place tube in incubator for 1-2 hours.
- 3. Mix 0.5 McFarland Standard and test suspension thoroughly with a vortexer or by shaking vigorously for 15 seconds.
- 4. In the presence of good lighting, visually compare turbidity of test suspension to the McFarland standard. This can be facilitated with the use of a Wickerham card or scale. (See illustration below.)
- 5. If test suspension appears too dense when compared to the McFarland standard, dilute the test suspension aseptically with additional broth or sterile saline.
- 6. If test suspension is too light, inoculate with additional harvested organisms or return to incubator for 1 hour then repeat procedure above.



## **Wickerham Turbidity Card**

If it is more difficult to see lines on Wickerham card through test suspension than through the McFarland standard, the suspension is too dense and must be diluted.

## (Printable Card)



Rev. 2/2014